

Quantitative Physical Stability Analysis of Paint, Ink and Coating Dispersions

Pascal Bru
Formulation

Formulation of Coatings

Introduction

- New paint, ink, and coating formulations rely on the formulator to provide optimal properties of the product throughout its entire lifetime.
- Pigments, binders, additives and solvents must be chosen and dispersed in such a way for the formulation to possess ideal properties including stability, viscosity, consistency, and overall quality.
- Many tools are available to chemists in order to analyze such qualities but some gaps remain. Extrapolation of results, correlation of measurements, and subjective analysis may improperly characterize the mixture.
- This can lead to poor performance of the product on an application side and the inability of a sample to redisperse after ageing, among other drawbacks.



Formulation of Coatings


Physical Stability Concerns

- A few issues are critical when discussing the physical stability and ageing/evolution of formulations related to inkjets, paints, and coatings.
 - Physical appearance (visual phase separation).
 - Redispersion capability once a sample has aged.
 - Prediction of sediment particle packing.
 - General shelf life.
- Additionally, R&D formulators would benefit from a fast method to **quantify** how changes in chemistry effect the formulation.
 - Regulations
 - Different raw material lots
 - New, improved formulations



Static Multiple Light Scattering Technology and Theory


Static Multiple Light Scattering can be used to monitor particle concentration and size changes inside of concentrated formulations.



Backscattering (BS) is a function of:

d : particle size

Φ : particle concentration



Repetition of the measurement provides:

Δd : change in particle size

$\Delta \Phi$: change in particle concentration

Allows for quantification and prediction of destabilization kinetics such as creaming, flocculation, and sedimentation.



SMLS Applications


Paint & Ink



Pigment Formulas

White Pigment Formulation Modification


- Three white pigment formulations were analyzed for top-layer clarification as an indicator of overall stability.
- This clarification kinetic can be quantified and the formulations ranked for stability well before any visual phase separation is detected by the naked eye.



Pigment Formulas

White Pigment Formulation Modification

- Light scattered in the clarification/phase separation zone can provide top-of-vial kinetic.
- This allows to quickly and quantitatively determine how the modification in the formula alters the final stability of product.




Pigment Formulas

White Pigment Formulation Modification

- Global stability kinetics (including any sedimentation and flocculation) show clear indication that formulation #3 exhibits greater instability than the other two.
- This formulation should be modified or the process changed to match the stability of the other two formulas.

Destabilisation Kinetics (Global)




Extremely useful data for formulation control

Inkjet Formulations

Effect of Pigment


- Inkjet samples are prone to destabilization and packing of the sediment. This hard packing can typically be redispersed by shaking or movement of a printing head.
- Kinetically, the packing is detected as a *decrease* in Δ BS. Even though there is an *increasing concentration*, the particles are *increasing in size* and dominating the signal change.



Inkjet Formulations

Effect of Pigment


- Samples with different pigments and formulations can be analyzed for the same types of destabilizations.
- This will allow for simple formulation adjustment and to see how different raw materials will behave in the formulation.



Inkjet Formulations

Effect of Pigment

- The mean value of light scattered at the bottom sediment layer is directly correlated to a packing kinetic.
- A few clicks will provide a detailed plot showing the progress of the packing. Note that the two samples are almost the same in this regard.



Advanced quantification of individual kinetics



Transparent Material Applications


Refinish Stability

- Clear coats, pre-treatments and topcoat formulations have a degree of transparency, allowing the particle motion of the system to be traced by the transmission detector of SMLS instrumentation.
- Three samples are tested at both RT (27 °C) and 50 °C to test for particle flocculation, sedimentation, and phase separation.



Transparent Material Applications

Refinish Stability – High T Acceleration



Flocculation is the first kinetic to be seen


Sedimentation causes a decrease in transmission

All samples exhibit similar behavior – but which are more stable?

Transparent Material Applications

Refinish Stability – High T Acceleration

- Bottom-of-sample analysis shows that flocculation dominates and is then taken over by sedimentation as the ΔT decreases.
- The slope of the graphs provides info into the sedimentation kinetic, even when other signals are interfering in the calculation zone.




Sample	Slope (18 hr to 3 days)
A, 50 °C (red)	-2.12%/day
B, 50 °C (blue)	-1.80%/day
C, 50 °C (green)	-9.00%/day

Transparent Material Applications

Refinish Stability – High T Acceleration

- Global stability kinetics show a distinct difference between samples analyzed at RT (but the naked eye cannot).
- Plots at 50 °C provide a clear distinction between the samples, which shows the same trend from the visual analysis of the samples.




3 days of analysis provide comparative stability data

Shelf Life Projections

Hard Settle Study: Toners

- Toners with the ability to hard pack over their shelf life were analyzed both with SMLS (5 days) as well as by standard methods (6 months).
- The results for both studies correlated exactly with the propensity of the formulations to not produce a hard pack, surely produce a hard pack, or to remain questionable for a hard pack prediction.




Decrease of BS at the bottom is evidence of packing

Large, broad clarifying layer (but still opaque)

Shelf Life Projections

Hard Settle Study: Toners

- In order to avoid the clarification kinetic in the global stability kinetic at the bottom of the sample (bottom 20% height) was analyzed.
- Hard packing is seen in the red (observed hard settle) and green (questionable hard settle) traces, evidenced by the high numbers.
- Samples that do not hard pack (blue trace) would have very low changes in this region.




6 month test completed in 5 days with the same results!

Cookware Coatings Applications

Stability of Sprayable Coatings

- White PTFE polymer dispersions for spray coatings show phase separations based upon the components added to the mixtures (analyzed at 45 °C).



Sample	Clarification intensity (Δ BS%, 36 hours)	Creaming phase volume (Δ BS%, 36 hours)	Total stability (36 hours)	Total stability (4 days)
1	-0.47	0.89	0.4	1.1
2	-0.37	0.70	0.3	0.6

Conclusions

Stability of Coating Applications

- Physical stability kinetics of formulations can be analyzed, quantified, and predicted using Static Multiple Light Scattering.
- Concentrated inkjet dispersions, paint and coating dispersions, and clear coat and refinish dispersion all benefit from this type of analysis.
- This allows for fast formulation optimization and accurate predictions into how the chemistry of a formulation affects its performance.

Thank you!
Please visit us at booth #61!

bru@formulation.com



www.formulation.com

STABILITY & SIZE

MICORHEOLOGY

RHEOLOGY ON CHIP